



ebioxVET



INFECTION CONTROL SYSTEM



ebioxVET is active worldwide and offers a range of products to veterinary practices where high level cleaning, decontamination and control of cross infection is required.

ebiox, advanced decontamination technology, have developed a comprehensive range of products specifically designed to remove the bio burden and eliminate micro organisms and harmful pathogens.

All products are manufactured to ISO 9002 and international quality standards and have been tested and approved.

Each product in the system has been individually designed for a specific purpose and represents the latest advance in cleaning and decontamination technology. When used separately each product will give outstanding results but when used together as part of a system the synergistic effect will minimise the potential for cross contamination.



Hospital Infection Research Laboratory, Birmingham (HIRL),

“ Thorough cleaning will remove micro organisms and the organic material on which they thrive. ebiox Ultra was the most effective product tested.”

ebioxVET INFECTION CONTROL SYSTEM

Ultra P

A patented solution to clean surgical instruments to molecular level*, removing tissue residues, fluids and other pathogenic debris. Use as a pre-soak, in ultrasonic baths or for manual cleaning of surgical instruments.

**Independent Test HIRL - removes more than 99.9% of pathogenic soiling.*

- Use to pre-soak instruments prior to autoclaving
- Particularly effective on difficult complex instruments
- Easily removes 'lubricating' jelly
- Safe to use on delicate devices, ophthalmic instruments and sensitive tubing
- Ideal for cold cleaning of endoscopes and heat sensitive instruments
- Pleasant fragrance, no release of irritating vapours



Presentation and Usage Directions

Super concentrated powder which quickly dissolves in warm water at a dilution rate of 20g per 2 litres of water. Open instruments to fullest extent and leave to soak for a period of 10 minutes. Remove and rinse

“ Conventional detergents and disinfectants simply do not have the power to deal with the rapidly changing bio-protective systems evolved by micro organisms.”

Order Code	Description	Case Quantity
EVUP6	Ultra P Surgical Instrument Cleaner - 600g	2 x 600g

Ultra SSC

Surgical Steel Reconditioning Powder - removes stains, 'rust' marks and dull shadow from the surface of surgical instruments and stainless steel. Ideal for reconditioning the damaged surface of surgical steel and reforming the bright protective metal film.

- Prolongs the investment in expensive surgical instruments
- Removes 'rust' and iron stains
- Removes other oxide films that cause dulling, turning instruments into "as new"

Presentation and Usage Directions

Supplied as a powder. Sprinkle onto a damp cloth and rub gently over the instrument and rinse.



Order Code	Description	Case Quantity
EVSC1	Ultra SSC Surgical Steel Reconditioning Powder - 100g	4 x 100g

Trionic

A convenient ready to use surface disinfectant and decontaminating spray based on Trionic advanced technology. Cleans, deodorises, disinfects and sanitises surfaces in one easy action. Designed for regular use on medical surfaces, cages, examination tables, trolleys and furniture.

- Convenient and effective
- Provides a clean and disinfected surface
- Total cell removal - no residual organisms left to develop into resistant strains
- Proven effective and fast acting
- Alcohol free - use on all washable surfaces including plastics
- Available in two presentations
- 5 litre presentation features 'calming' fragrance
- Spray presentation features 'medical' fragrance

Presentation and Usage Directions

Apply with a spray, mop or cloth onto any water safe surface. Leave for a few minutes, then wipe, mop or rinse thoroughly.



Order Code	Description	Case Quantity
EVTP5	Trionic Concentrate - 5 litre	2 x 5 litre
EVTH5	Trionic Spray Disinfectant - 500ml	12 x 500ml

Trionic Plus

Cleaning

Trionic Plus is a new technology surface agent which will penetrate and completely remove the biofilm and embedded pathogens. Trionic Plus is probably the safest and most powerful biocide available.

Unique patented formula to clean, deodorise, disinfect and sanitise surfaces and maintain a hygienic environment.

- **Kills Parvo Virus in seconds**
- Immediately eliminates bacteria, viruses, fungi and spores
- More effective than aldehydes, acids, chlorine and alcohol
- User and environmentally friendly
- Use with confidence in kennels, consulting rooms, floors, walls, cages, tables, furniture and public areas
- Uniquely cleans, deodorises, disinfects and sanitises in one action
- Special 'calming' fragrance

Hospital Infection Research Laboratory, Birmingham,

“**Trionic Plus passed EN1276 achieving greater than 99.999% reduction of all test organisms within 30 seconds.**”

Presentation and Usage Directions

Trionic Plus is a two part solution that requires activation before use. To prepare approximately 1 litre of active disinfectant solution - add 2 capfuls (2 x 30ml) of concentrate Trionic Plus to a mixing vessel. Add 1 capful (1 x 30ml) of activator to the same mixing vessel. Gently stir the two liquids together and leave for exactly 5 minutes. Add this concentrate solution to 1 litre of water. To prepare greater volumes of disinfectant follow the same procedure using volumes indicated in the following table.

Volume of Disinfectant	No of Caps of Trionic Plus	No of Caps of Trionic Activator	Volume of Water
2.18 litre	4	2	2 litre
6.54 litre	12	6	6 litre
10.9 litre	20	10	10 litre

Apply with a spray, cloth, sluice or mop onto any water safe surface, then wipe clean and rinse surface. For maximum disinfection leave on surface for up to 10 minutes, then importantly, thoroughly rinse or flush surface with water.



Order Code	Description	Case Quantity
EVTK15	Trionic Plus and Activator Kit	4 x 1.5 litre Plus 2 x 1.5 litre Activator
EVTP15	Trionic Plus	6 x 1.5 litre
EVTA15	Activator	6 x 1.5 litre

Dr Paul Burr MRCVS - Veterinary Virologist,
Pentlands Science Park, Midlothian,

“**Trionic Plus has a significant effect on Canine Parvo Virus. The effect seen was rapid, within one minute no virus was detected.**”



Esense Hand Scrub

Effective hand washing is one of the most vital infection control measures. Esense Hand Scrub from Ebiox is a highly effective anti-bacterial cleanser which contains moisturisers and conditioners to protect the hands.

Highly efficient in removing the transient organisms as well as reducing resident organisms to a 'safe' level. Effective against gram-negative bacteria such as Staphylococcus aureus, and E-Coli 157; antibiotic resistant bacteria such as MRSA and VRE; medically important viruses including Rotavirus, Vaccinia, Flu A2 Hong Kong and Herpes simplex type.

- Moisturising leaving hands feeling soft and conditioned
- Chlorhexadine and alcohol free to avoid skin irritation and dryness
- Highly effective anti microbial protection
- Convenient pump dispenser
- CEN Approved passing the EN1499 and prEN12054 tests

Presentation and Usage Directions

Apply recommended 3ml quantity to hands and wash well with water using approved technique (ie Aycliffe technique in palm to palm, left palm over right dorsum and reverse, fingers interlock, back of fingers, thumbs and wrists). Rinse thoroughly and dry hands in suitable manner.

Use regularly throughout the day.

Order Code	Description	Quantity
EVHS6	Esense Hand Scrub - 600ml	4 x 600ml bottles with 2 pump dispensers (2 x 600ml with 2 refills)

Hospital Infection Research Laboratory, Birmingham,

“**Esense Hand Scrub from ebiox conforms to European efficacy standard EN1499 and is significantly superior to the reference product.**”



Esense Hand Disinfectant

Scientifically advanced hand lotion disinfectant to assist in the control of cross infection from hand contact.

Mild lotion with low skin irritancy, alcohol free and moisturising essential for frequent use applications.

Effective against a wide range of both gram-positive and gram-negative bacteria such as Staphylococcus aureus, and E-Coli 157; antibiotic resistant bacteria such as MRSA and VRE; medically important viruses including Rotavirus, Vaccinia, Flu A2 Hong Kong and Herpes simplex type.

- Moisturising leaving hands feeling soft and conditioned
- Highly effective anti microbial protection
- Convenient pump dispenser
- CEN Approved achieving greater than 5 log reduction in skin micro organisms and passing prEN12054 test

Presentation and Usage Directions

Apply recommended 3ml quantity to hands. Apply using Aycliffe technique in palm to palm, left palm over right dorsum and reverse, fingers interlock, back of fingers, thumbs and wrists.

Use regularly throughout the day.

Order Code	Description	Case Quantity
EVHD6	Esense Hand Disinfectant - 600ml	4 x 600ml bottles with 2 pump dispensers (2 x 600ml with 2 refills)

Hospital Infection Research Laboratory, Birmingham,

“**Esense Hand Disinfectant is a moisturising lotion which is proven as effective as alcohol but without causing skin dryness or irritation.**”



ebioxVET

INFECTION CONTROL SYSTEM

ebioxVET

1 Whitehall

Whitehall Road

Leeds

LS1 4HR

T: +44(0)113 390 6060

F: +44(0)113 390 6100

E: sales@ebioxvet.co.uk

W: www.ebioxvet.co.uk

ebioxVET is a division of Ebiox Limited



Advanced Decontamination Technology